# **ResQProducts**

SPECIALISING IN MOBILITY IMPAIRED EVACUATION

### **PRODUCT DATA SHEET**

## ResQmat

#### DESCRIPTION

The ResQMat enables carers to rapidly and safely evacuate people with mobility problems in an emergency. The built in padding and extra pillow under the head mean that the ResQmat is safe and comfortable for people when going down stairs and over rough surfaces.

With three wide Velcro straps patients being moved are completely enclosed, meaning the patient is contained once they are in a safe location, which is especially important when working with people with dementia who may wander.

#### **FEATURES**

- 300kg Weight capacity
- 2 person operation
- Can be moved easily over virtually any surface, especially over rough surfaces
- Narrow width & flexibility
- Wipe clean storage bag
- Can be used down stairs



#### SPECIFICATION

Height	000mm (folded)
Width	000mm (folded)
Depth	000mm (folded)
Weight	00kg (unladen)
Payload Capacity	300kg

(EVAC+CHAIR

INTERNATIONAL

The ResQmat is made under licence from Pelican Manufacturing P/L, in the U.K. by Evac+Chair International Ltd.

Evac+Chair International, ParAid House, Weston Lane, Birmingham, B11 3RS, United Kingdom T: +44 (0)121 706 6744 **W:** www.evacchair.co.uk

